

Take the Dako Challenge

You already use our kappa and lambda light chain conjugates.
Now, sample another reagent in your flow cytometry lab.
You'll be back for more.



With more than 40 years experience in the diagnostics industry specializing in antibody development and manufacturing, Dako has a worldwide reputation for creating products of first-class functionality, quality and reliability.

This expertise can be found in our comprehensive line of more than 240 reagents for flow cytometry including our world-class polyclonal kappa and lambda light chain conjugates.

As your preferred supplier of flow products, we guarantee:

- Prompt delivery to ensure uninterrupted productivity
- A technical support team dedicated to keeping your laboratory up and running
- Premium high-quality flow products with lot-to-lot consistency

Be sure to ask for an evaluation sample of a high-quality Dako reagent for your flow cytometry lab today. To find out more about this limited time offer or our comprehensive product portfolio, please contact your local distributor.

Code	Source	Product
------	--------	---------

Kappa Light Chain Conjugates

C0222	Rb a Hu	Kappa Light Chains/APC, Rabbit F(ab') ₂ ★, for Flow Cytometry
F0434	Rb a Hu	Kappa Light Chains/FITC, Rabbit F(ab') ₂ ★, for Flow Cytometry
R0436	Rb a Hu	Kappa Light Chains/RPE, Rabbit F(ab') ₂ ★, for Flow Cytometry
FR048	a Hu	Kappa Light Chains/FITC, Rabbit F(ab') ₂ ★ + CD19/RPE, Clone HD37, MultiMix™ Dual-Color, for Flow Cytometry
FR481	Rb a Hu	Kappa Light Chains/FITC, Rabbit F(ab') ₂ ★ + Lambda Light Chains/RPE, Rabbit F(ab') ₂ ★, MultiMix™ Dual-Color, for Flow Cytometry
TC051	a Hu	Kappa Light Chains/FITC, Rabbit F(ab') ₂ ★ + Lambda Light Chains/RPE, Rabbit F(ab') ₂ ★ + CD19/RPE-Cy5, Clone HD37, MultiMix™ Triple-Color, for Flow Cytometry
TC669	a Hu	Kappa Light Chains/APC, Rabbit F(ab') ₂ ★ + CD19/FITC, Clone HD37 + Lambda Light Chains/RPE, Rabbit F(ab') ₂ ★, MultiMix™ Triple-Color, for Flow Cytometry
TC670	a Hu	Kappa Light Chains/APC, Rabbit F(ab') ₂ ★ + Lambda Light Chains/RPE, Rabbit F(ab') ₂ ★ + Plasma Cell/FITC, Clone VS38c, MultiMix™ Triple-Color, for Flow Cytometry

Lambda Light Chain Conjugates

F0435	Rb a Hu	Lambda Light Chains/FITC, Rabbit F(ab') ₂ ★, for Flow Cytometry
R0437	Rb a Hu	Lambda Light Chains/RPE, Rabbit F(ab') ₂ ★, for Flow Cytometry
FR044	a Hu	Lambda Light Chains/FITC, Rabbit F(ab') ₂ ★ + CD19/RPE, Clone HD37, MultiMix™ Dual-Color, for Flow Cytometry
FR481	Rb a Hu	Lambda Light Chains/RPE, Rabbit F(ab') ₂ ★ + Kappa Light Chains/FITC, Rabbit F(ab') ₂ ★, MultiMix™ Dual-Color, for Flow Cytometry
TC669	a Hu	Lambda Light Chains/RPE, Rabbit F(ab') ₂ ★ + CD19/FITC, Clone HD37 + Kappa Light Chains/APC, Rabbit F(ab') ₂ ★, MultiMix™ Triple-Color, for Flow Cytometry
TC051	a Hu	Lambda Light Chains/RPE, Rabbit F(ab') ₂ ★ + Kappa Light Chains/FITC, Rabbit F(ab') ₂ ★ + CD19/RPE-Cy5, Clone HD37, MultiMix™ Triple-Color, for Flow Cytometry
TC670	a Hu	Lambda Light Chains/RPE, Rabbit F(ab') ₂ ★ + Kappa Light Chains/APC, Rabbit F(ab') ₂ ★ + Plasma Cell/FITC, Clone VS38c, MultiMix™ Triple-Color, for Flow Cytometry

* Affinity - isolated antibody



Corporate Headquarters
Denmark
+45 44 85 95 00

Find your local distributor through:

www.dako.com

or contact
Reagent Partnership Division on
reagent.partnerships@dako.com

Distributors in more
than 50 countries